

# 3-PLY MEDICAL FACE MASKS

## TYPE IIR BFE ≥ 98%



**With your own logo and coloured elastics**



**Product:**

- ✓ Medical face mask
- ✓ 17.5 cm x 9.5 cm
- ✓ Medical device of class I in application of annex IX , European directive 93/42/CEE
- ✓ Typ IIR according to EN14683 + AC 2019, CE certified
- ✓ Mask can be marked with a logo

**Technical data sheet:**

|  | Norm specification      | Mask performance      |
|--|-------------------------|-----------------------|
| <b>Bacterial filtration efficiency</b> | ≥ 98%                   | ≥ 99.5%               |
| <b>Respirability</b>                   | < 60 Pa/cm <sup>2</sup> | 48 Pa/cm <sup>2</sup> |
| <b>Microbial cleanliness</b>           | < 30 CFU/g              | < 30 CFU/g            |
| <b>Splash resistance</b>               | >16 kPa                 | >16 kPa               |

Logo to be supplied under following formats : Bmp / jpg ; jif ; jpeg / pcx / png / tif ; tiff / tga.

**Packaging:**

**1. Box of 50 masks:**



72 boxes (3.600 masks) are supplied in cardboard boxes (575 x 420 x 620 mm).  
 8 cardboard boxes as a maximum per pallet (28.800 masks).  
 Min quantity to order 3.600 masks.

**2. Plastic bag of 10 masks**

**Glanzstoff**  
**TEXTILCORD**  
An Indorama Ventures Company

**MD Medical face masks**

- EN 14683 + AC 2019
- Type IIR BFE >=98%
- Professional
- High breathable
- Comfortable
- High elasticity ear loop

**We protect you**

**10 pcs**

17.5x9.5cm

330 Plastics (3.300 masks) are supplied in cardboard boxes (575 x 420 x 620 mm).  
 8 cardboard boxes as a maximum per pallet (26.400 masks).  
 Min quantity to order 3.300 masks.

# 3-PLY MEDICAL FACE MASKS

## TYPE IIR BFE ≥ 98%



## PURCHASE ORDER

| NAME              |  |
|-------------------|--|
| Company Name      |  |
| Address           |  |
| Country           |  |
| CONTACT           |  |
| Phone Number      |  |
| Fax Number        |  |
| Contact           |  |
| E-mail            |  |
| ORDER INFORMATION |  |

Date: \_\_\_\_\_

Expected delivery date: \_\_\_\_\_

| Masks quantity | Internal code | 50 masks / boxes | 10 masks / plastic bag | With logo | unit price/mask without TVA<br>FCA Steinfort | Mask quantity ordered | Total price |
|----------------|---------------|------------------|------------------------|-----------|--|-----------------------|-------------|
| 3.600          | PCM0003E      | yes              | no                     | yes       | 0.12€  |                       |             |
| 3.600          | PCM0003B      | yes              | no                     | no        | 0.10€  |                       |             |
| 3.300          | PCM0003F      | no               | yes                    | yes       | 0.14€  |                       |             |
| 3.300          | PCM0003C      | no               | yes                    | no        | 0.13€  |                       |             |

Please send your Purchase Order to: [mask.textilcord@glanzstoff.com](mailto:mask.textilcord@glanzstoff.com)