

3-PLY MEDICAL FACE MASKS

TYPE IIR BFE ≥ 98%



Storage conditions:



Product validity: 2 years

Instruction for use:



1. Open the mask



2. Nose strip at the top, hung on the ear



3. Adjust the mask size to completely cover the nose and the mouth



4. Use both hands to regulate the binding bridge on the nose on both sides of the masks



5. Stamped side of mask to be placed outside

This product is a single-use product.

It is strictly forbidden to clean or reuse it.

Disposable face mask to be replaced every 4 hours.